

Data Sheet

HUMAN NATRIURETIC PEPTIDE RECEPTOR B [NPR-B] (aa 288-302)

SYNTHETIC

Catalog No.:	AE1032AG.1 / AE1032AG.2
Name:	Human NPR-B (aa 288-302)
Synonyms:	Atrial natriuretic peptide receptor 2, ANP-B, ANPR-B, Guanylate cyclase B (GC-B), EC=4.6.1.2
Sequence:	QNRLLRAREDFGVE
Swiss-Prot No:	P20594
Gene Information:	Gene Name: NPR2, ANPRB GeneID: 4882
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	18 May 2022

For research use only

Publishing research using AE1032AG? Please let us know so that we can cite your publication as a reference.

