

Data Sheet

HUMAN NATRIURETIC PEPTIDE RECEPTOR C [NPR-C] (aa 199-213)

SYNTHETIC

Catalog No.:	AE1030AG.1 / AE1030AG.2
Name:	Human NPR-C (199-213) (QEGLHTSIYSFDET)
Synonyms:	Atrial natriuretic peptide receptor 3, Atrial natriuretic peptide clearance receptor
Sequence:	QEGLHTSIYSFDET
MW:	1.65 kDa
Swiss-Prot No:	P17342
Gene Information:	Gene Name: NPR3, ANPRC GeneID: 4883
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20°C Repeated thawing and freezing must be avoided
Last updated on:	12 May 2022

For research use only

Publishing research using AE1030AG? Please let us know so that we can cite your publication as a reference.

