

Data Sheet

MOUSE NIDOGEN 1

RECOMBINANT

Catalog no.:	AP1003AG.1 / AP1003AG.2
Name:	Mouse Nidogen 1
Synonyms:	Entactin
Swiss-Prot No:	P10493
Gene Information:	Gene Name: Nid1 GeneID: 18073
Purity:	> 95% (SDS-PAGE)
Matrix:	50 mM HEPES, 150 mM NaCl, pH 7.5
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	20 May 2022

For research use only

Publishing research using AP1003AG? Please let us know so that we can cite your publication as a reference.