



Data Sheet

NORMAL GUINEA PIG SERUM

Catalog no.:	NS4003
Host:	Guinea Pig
Matrix:	Serum, 137 mM D-(-)-mannitol
Contents:	5 ml (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	26 March 2025

For research use only

Publishing research using NS4003? Please let us know so that we can cite your publication as a reference.

