



Data Sheet

NORMAL GUINEA PIG IgG

Catalog no.:	NG4003.1 / NG4003.2
Host:	Guinea Pig
Matrix:	IgG, Protein A purified from serum; PBS pH 7.2, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol
Contents:	1 mg / 5 mg (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	2 - 8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
Last updated on:	27 March 2025

For research use only

Publishing research using NG4003? Please let us know so that we can cite your publication as a reference.

