



## Data Sheet

# NORMAL CHICKEN IgY

<b>Catalog no.:</b>	NG4011.1 / NG4011.2
<b>Host:</b>	Chicken
<b>Matrix:</b>	IgY, purified from egg yolk; PBS pH 7.2
<b>Contents:</b>	1 mg / 5 mg (lyophilized) Resuspend in aqua bidest for required concentration.
<b>Store at:</b>	2 - 8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
<b>Last updated on:</b>	4 May 2022

### For research use only

**Publishing research using NG4011? Please let us know so that we can cite your publication as a reference.**

