



Data Sheet

NORMAL RABBIT IgG

Catalog no.:	NG4004.1 / NG4004.2
Host:	Rabbit
Matrix:	IgG, Protein A purified from serum; PBS pH 7.2
Contents:	1 mg / 5 mg (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	2 - 8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
Last updated on:	4 July 2024

For research use only

Publishing research using NG4004? Please let us know so that we can cite your publication as a reference.

