



Data Sheet

NORMAL CHICKEN IgY

Catalog no.:	NG4011.1 / NG4011.2
Host:	Chicken
Matrix:	IgY, Protein G purified from egg yolk; PBS pH 7.2, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol
Contents:	1 mg / 5 mg (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	2 - 8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
Last updated on:	26 March 2025

For research use only

Publishing research using NG4011? Please let us know so that we can cite your publication as a reference.

