



Data Sheet

NORMAL GOAT IgG

Catalog no.:	NG4002.1 / NG4002.2
Host:	Goat
Matrix:	IgG, Protein G purified from serum; PBS pH 7.2
Contents:	1 mg / 5 mg (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	2 - 8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
Last updated on:	4 May 2022

For research use only

Publishing research using NG4002? Please let us know so that we can cite your publication as a reference.

