



## Data Sheet

### NORMAL RABBIT IgG

<b>Catalog no.:</b>	NG4004.1 / NG4004.2
<b>Host:</b>	Rabbit
<b>Matrix:</b>	IgG, Protein G purified from serum; PBS pH 7.2, 0.08% NaN <sub>3</sub>
<b>Contents:</b>	1 mg / 5 mg (lyophilized) Resuspend in aqua bidest for required concentration.
<b>Store at:</b>	2 - 8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
<b>Last updated on:</b>	4 May 2022

**For research use only**

**Publishing research using NG4004? Please let us know so that we can cite your publication as a reference.**

