



Data Sheet

HUMAN OSTEOINDUCTIVE FACTOR [OIF] (aa 107-121)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1051.1 / AE1051.2
Immunogen:	Synthetic human OIF (aa 107-121) bTG-conjugated (RFNKIKKLTAKDFAD)
Synonyms:	Mimecan, Osteoglycin
Swiss-Prot No:	P20774
Gene Information:	Gene Name: OGN, OIF, SLRR3A GenelD: 4969
Host:	Rabbit
Matrix:	Serum
Specificity:	Human OIF (aa 107-121) There was no cross-reactivity obtained with human OIF (aa 187-198), human Osteopontin, and human Bone Sialoprotein.
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE1051? Please let us know so that we can cite your publication as a reference.

