



Data Sheet

RAT OSTEOACTIVIN (aa 551-568)

ANTIBODY, POLYCLONAL

Catalog no.:	A 1080.1 / A 1080.2
Immunogen:	Synthetic rat osteoactivin (aa 551-568) (KAPFSRGDREKDPLLQDK)
Swiss-Prot No:	Q9QZF6
Gene Information:	Gene Name: Gpnmb GenelD: 113955
Host:	Rabbit
Matrix:	Serum
Specificity:	Rat osteoactivin (aa 551-568)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1080? Please let us know so that we can cite your publication as a reference.

