

## Data Sheet

# HUMAN OSTEOCALCIN (aa 1-22)

## SYNTHETIC

<b>Catalog no.:</b>	A 4211AG.1 / A 4211AG.2
<b>Name:</b>	Human Osteocalcin (aa 1-22) (YLYQWLGAPVPYPDPLEPRREV)
<b>Synonyms:</b>	Bone Gla protein (BGP), Gamma-carboxyglutamic acid-containing protein
<b>Swiss-Prot No:</b>	P02818
<b>Gene Information:</b>	Gene Name: BGLAP GenelD: 632
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	04 May 2022

**For research use only**

**Publishing research using A 4211AG? Please let us know so that we can cite your publication as a reference.**

