



Data Sheet

HUMAN OSTEONECTIN (aa 272-286)

SYNTHETIC

Catalog no.:	A 4233AG.1 / A 4233AG.2
Name:	Human Osteonectin (aa 272-286)
Synonyms:	Secreted protein acidic and rich in cysteine (SPARC), Basement-membrane protein 40 (BM-40)
Sequence:	CFGIKQKDIDKDLVI
MW:	1.74 kDa
Swiss-Prot No:	P09486
Gene Information	Gene name: SPARC, ON GenelD: 6678
Source:	Synthetic
Purity:	> 95% (SDS-PAGE)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 0.1 % acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	5 May 2022

For research use only

Publishing research using A 4233AG? Please let us know so that we can cite your publication as a reference.

