

Data Sheet

HUMAN OSTEONECTIN

RECOMBINANT

Catalog no.:	A 4225AG.1 / A 4225AG.2
Name:	Human Osteonectin
Synonyms:	SPARC
Swiss-Prot No:	P09486
Gene Information	Gene name: SPARC GeneID: 6678
Source:	Produced serum-free in a human cell line
Purification:	Affinity chromatography
Storage buffer:	PBS pH 7.4
Purity:	> 95% (SDS-PAGE)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Reference:	1. Salza R, Peysselon F, Chautard E, Faye C, Moschovich L, Weiss T, Perrin-Cocon L, Lotteau V, Kessler E, Ricard-Blum S (2014) Extended interaction network of procollagen C-proteinase enhancer-1 in the extracellular matrix. <i>Biochem J</i> 457(1): 137-149
Last updated on:	05 May 2022

For research use only

Publishing research using A 4225AG? Please let us know so that we can cite your publication as a reference.

