



Data Sheet

HUMAN OSTEOPONTIN (aa 59-74)

SYNTHETIC

Catalog no.:	A 4226AG.1 / A 4226AG.2
Name:	Human Osteopontin (aa 59-74)
Synonyms:	Bone sialoprotein 1, Secreted phosphoprotein 1(SPP-1), Urinary stone protein, Nephropontin, Uropontin
Sequence:	PSKSNESHDMDDMDD
MW:	1.86 kDa
Swiss-Prot No:	P10451
Gene Information:	Gene Name: SPP1, BNSP, OPN GenelD: 6696
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	05 May 2022

For research use only

Publishing research using A 4226AG? Please let us know so that we can cite your publication as a reference.

