

## Data Sheet

# HUMAN OSTEOPONTIN (aa 164-179)

## SYNTHETIC

<b>Catalog no.:</b>	A 4227AG.1 / A 4227AG.2
<b>Name:</b>	Human Osteopontin (aa 59-74)
<b>Synonyms:</b>	Bone sialoprotein 1, Secreted phosphoprotein 1 (SPP-1), Urinary stone protein, Nephropontin, Uropontin
<b>Sequence:</b>	QYPDATDEDITSHMES
<b>MW:</b>	1.84 kDa
<b>Swiss-Prot No:</b>	P10451
<b>Gene Information:</b>	Gene Name: SPP1, BNSP, OPN GeneID: 6696
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	5 May 2022

**For research use only**

**Publishing research using A 4227AG? Please let us know so that we can cite your publication as a reference.**

