



Data Sheet

CHICKEN OVALBUMIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 1222.1 / A 1222.2
Immunogen:	Ovalbumin purified from chicken egg whites
Host:	Goat
Matrix:	Serum
Specificity:	Chicken ovalbumin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: ELISA (1:500)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); -20 °C (dissolved)
Repeated thawing and freezing must be avoided

Last updated on: 16 March 2022

For research use only

Publishing research using A 1222? Please let us know so that we can cite your publication as a reference.

