



Data Sheet

HUMAN PANCREATIC ELASTASE-3B

ANTIBODY, MONOCLONAL

Catalog no.:	ADR1001.1/ADR1001.2
Immunogen:	Recombinant human Elastase-3B
Synonyms:	CELA 3B, Chymotrypsin-like elastase family member 3B, EC:3.4.21.70, Elastase IIIB, Protease E
Swiss-Prot No:	P08861
Gene Information:	Gene Name: CELA3B, ELA3B Gene ID: 23436
Host:	Mouse Balb/c
Clone no.:	FP112
Isotype:	IgG1 kappa
Matrix:	PBS pH 7.4, 0.09% NaN ₃
Specificity:	Human fecal Elastase
Contents:	10 µg / 100 µg (liquid)
Known applications:	ELISA (coating antibody)
	<small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C
Last updated on:	14 June 2024

For research use only

Publishing research using ADR1001? Please let us know so that we can cite your publication as a reference.

