



Data Sheet

HUMAN PANCREATIC ELASTASE-3B RECOMBINANT

Catalog no.:	ADR1005AG.1 / ADR1005AG.2
Name:	Human Elastase-3B
Synonyms:	Chymotrypsin-like elastase family member 3B, EC:3.4.21.70, Elastase IIIB, Protease E
MW:	28.4 kDa
Swiss-Prot No:	P08861
Gene Information:	Gene Name: CELA3B, ELA3B Gene ID: 23436
Source:	Mammalian cell line
Purification:	Immobilized metal ion affinity chromatography (IMAC), preparative SEC
Storage buffer:	PBS pH 7.4
Purity:	>90%
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	12 May 2022

For research use only

Publishing research using ADR1005AG? Please let us know so that we can cite your publication as a reference.

