



Data Sheet

HUMAN PARATHYROID HORMONE (PTH) (aa 1-10)

ANTIBODY, MONOCLONAL

Catalog no.:	AK1109.1 / AK1109.2
Immunogen:	Synthetic human PTH peptide (aa 1-10) poly-Lysine conjugated
Synonyms:	Parathyrin, Parathormone
Swiss-Prot No:	P01270
Gene Information:	Gene name: PTH GeneID: 5741
Host:	Mouse Balb/c
Clone no.:	L7-12A-F12
Isotype:	IgG1 _{kappa}
Matrix:	Cell supernatant, Protein G purified; PBS pH 7.4
Specificity:	Human PTH peptide <1p mol/l (aa 1-3; 1-10; 1-34; 1-37; 1-84)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA (50 ng/ml); ILMA (50 ng/ml) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	27 April 2022

For research use only

Publishing research using AK1109? Please let us know so that we can cite your publication as a reference.

