



Data Sheet

HUMAN PARATHYROID HORMONE [PTH]

(aa 53-84)

ANTIBODY, POLYCLONAL

Catalog no.:	A1113.1 / A1113.2
Immunogen:	Synthetic human PTH (aa 53-84) bTG- conjugated (KKEDNVLVESHEKSLGEADKADVNVLTKAKSQ)
Synonyms:	Parathyrin, Parathormone
Swiss-Prot No:	P01270
Gene Information:	Gene name: PTH GenelD: 5741
Host:	Goat
Matrix:	Serum
Specificity:	Human PTH (aa 53-68; 53-84; 1-84)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:10 000); ELISA (1:5 000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1113? Please let us know so that we can cite your publication as a reference.

