



Data Sheet

RAT PARATHYROID HORMONE (PTH) (aa 1-34)

ANTIBODY, POLYCLONAL

Catalog no.: A 1711.1 / A 1711.2

Immunogen: Synthetic rat PTH (aa 1-34) bTG- conjugated
(AVSEIQLMHNLGKHLASVERMQWLRKKLQDVHNF)

Synonyms: Parathyrin, Parathormone

Swiss-Prot No: P01270

Gene Information: Gene Name: PTH
GeneID: 5741

Host: Goat

Matrix: Serum

Specificity: Rat PTH (aa 1-34)

There was no cross reactivity obtained with rat PTH (aa 53-84).

Contents: 20 µl / 100 µl (lyophilized)

Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: RIA (1:2000)¹

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

References: 1. Klöppinger, M., Armbruster, F. P., Schaefer, L., Tampe, J., and Schmidt-Gayk, H. Homologer N-terminaler Radioimmunoassay für Ratten-Parathormon(1-34). In: Atkinson, M., Raue, F., Schmidt-Gayk, H., Schweickert, H. U., Wuster, C., and Ziegler, R. 22. 1992. Osteologie '92: Jahrestagung der Sektion Calcium-regulierende Hormone und Knochenstoffwechsel (CRHUKS) der Deutschen Gesellschaft für Endokrinologie, 1-2 Oktober 1992, Heidelberg. 1-1-0092.

Last updated on: 16 March 2022

For research use only

Publishing research using A 1711? Please let us know so that we can cite your publication as a reference.

