



## Data Sheet

# HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 1-34)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 7881.3 / A 7881.4
<b>Immunogen:</b>	Synthetic human PTHrP (aa 1-34) KLH conjugated (AVSEHQLLHDKGKSIQDLRRRFFLHHLIAEIHTA)
<b>Synonyms:</b>	Parathyroid hormone-like protein (PLP)
<b>Swiss-Prot No:</b>	P12272
<b>Gene Information:</b>	Gene Name: PTHLH, PTHRP GenelD: 5744
<b>Host:</b>	Rabbit
<b>Matrix:</b>	IgG, Protein A purified from serum; PBS pH 7.4, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol
<b>Specificity:</b>	Human PTHrP (aa 1-34)
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	13 March 2025

### For research use only

**Publishing research using A 7881? Please let us know so that we can cite your publication as a reference.**

