



Data Sheet

HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 53-86)

ANTIBODY, POLYCLONAL

Catalog no.:	A 1118.1 / A 1118.2
Immunogen:	Synthetic human PTHrP (aa 53-86) KLH-conjugated (KNHPVRFSGSDDEGRYLTQETNKVETYKEQPLKTP)
Synonyms:	Parathyroid hormone-like protein (PLP)
Swiss-Prot No:	P12272
Gene Information:	Gene Name: PTHLH, PTHRP GenelD: 5744
Host:	Goat
Matrix:	Serum
Specificity:	Human PTHrP (aa 53-86; 1-86)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:2000) ¹ This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Blind E, Raue F, Gotzmann J, Schmidt-Gayk H, Kohl B, Ziegler R (1992). Circulating levels of midregional parathyroid hormone-related protein in hypercalcaemia of malignancy. <i>Clin Endocrinol (Oxf)</i> 37 (3): 290-297.
Last updated on:	16 March 2022

For research use only

Publishing research using A 1118? Please let us know so that we can cite your publication as a reference.

