



## Data Sheet

# HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 1-34)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 4243.1 / A 4243.2
<b>Immunogen:</b>	Synthetic human PTHrP (aa 1-34) KLH conjugated (AVSEHQLLHDKGKSIQDLRRRFFLHHLIAEIHTA)
<b>Synonyms:</b>	Parathyroid hormone-like protein (PLP)
<b>Swiss-Prot No:</b>	P12272
<b>Gene Information:</b>	Gene Name: PTHLH, PTHRP GenelD: 5744
<b>Host:</b>	Goat
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human PTHrP (aa 1-34)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA (1:2000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	6 April 2022

### For research use only

**Publishing research using A 4243? Please let us know so that we can cite your publication as a reference.**

