



Data Sheet

HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 1-34)

ANTIBODY, POLYCLONAL

Catalog no.:	A 7881.3 / A 7881.4
Immunogen:	Synthetic human PTHrP (aa 1-34) KLH conjugated (AVSEHQLLHDKGKSIQDLRRRFFLHHLIAEIHTA)
Synonyms:	Parathyroid hormone-like protein (PLP)
Swiss-Prot No:	P12272
Gene Information:	Gene Name: PTHLH, PTHRP GenelD: 5744
Host:	Rabbit
Matrix:	IgG, Protein G purified from serum; PBS pH 7.4
Specificity:	Human PTHrP (aa 1-34)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	RIA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	8 April 2022

For research use only

Publishing research using A 7881? Please let us know so that we can cite your publication as a reference.

