



Data Sheet

HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 1-34)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | A 4243.3 / A 4243.4 |
| Immunogen: | Synthetic human PTHrP (aa 1-34) KLH conjugated (AVSEHQLLHDKGKSIQDLRRRFFLHHLIAEIHTA) |
| Synonyms: | Parathyroid hormone-like protein (PLP) |
| Swiss-Prot No: | P12272 |
| Gene Information: | Gene Name: PTHLH, PTHRP GenelD: 5744 |
| Host: | Goat |
| Matrix: | IgG, Protein G purified from serum; PBS pH 7.4, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol |
| Specificity: | Human PTHrP (aa 1-34) |
| Contents: | 10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest. |
| Known applications: | RIA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 13 March 2025 |

For research use only

Publishing research using A 4243? Please let us know so that we can cite your publication as a reference.

