

## Data Sheet

# PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 107-111)

## SYNTHETIC

<b>Catalog no.:</b>	A 4241AG.1 / A 4241AG.2
<b>Name:</b>	Human PTHrP (107-111)
<b>Synonyms:</b>	Parathyroid hormone-like protein (PLP), Osteostatin
<b>Sequence:</b>	TRSAW
<b>MW:</b>	0.55 kDa
<b>Swiss-Prot No:</b>	P12272
<b>Gene Information</b>	Gene name: PTHLH, PTHRP GeneID: 5744
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	5 May 2022

### For research use only

**Publishing research using A 4241AG? Please let us know so that we can cite your publication as a reference.**

