

Data Sheet

HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 1-34)

SYNTHETIC

Catalog No.:	A 7881AG.1 / A 7881AG.2
Name:	Human PTHrP (aa 1-34)
Synonyms:	Parathyroid hormone-like protein (PLP)
Sequence:	AVSEHQLLHDKGKSIQDLRRRFFLHHLIAEIHTA
MW:	4.02 kDa
Swiss-Prot No:	P12272
Gene Information:	Gene Name: PTHLH, PTHRP GeneID: 5744
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	6 May 2022

For research use only

Publishing research using A 7881AG? Please let us know so that we can cite your publication as a reference.

