



Data Sheet

HUMAN PMN-ELASTASE

ANTIBODY, MONOCLONAL

Catalog no.:	A 6842.1 / A 6842.2
Immunogen:	Human PMN-Elastase
Synonyms:	Neutrophil elastase, Bone marrow serine protease, Elastase-2, Human leukocyte elastase (HLE), Medullasin, EC:3.4.21.37
Swiss-Prot No:	P08246
Gene Information:	Gene Name: ELANE, ELA2 GeneID: 1991
Host:	Mouse Balb/c
Clone no.:	IDmab-001
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Protein G purified, PBS pH 7.4, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol
Specificity:	Free/non-complexed PMN-Elastase, PMN-Elastase complexed to alpha1 proteinase inhibitor
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8°C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	29 May 2024

For research use only

Publishing research using A 6842? Please let us know so that we can cite your publication as a reference.

