



Data Sheet

HUMAN PMN ELASTASE

ANTIBODY, POLYCLONAL

Catalog no.:	A 6841.1 / A 6841.2
Immunogen:	Human purified PMN Elastase
Synonyms:	Neutrophil elastase, Bone marrow serine protease, Elastase-2, Human leukocyte elastase (HLE), Medullasin
Swiss-Prot No:	P08246
Gene Information:	Gene Name: ELANE, ELA2 GeneID: 1991
Host:	Goat
Matrix:	Serum
Specificity:	Human PMN Elastase
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:8000), RIA (1: 10 000), Western Blot (1:500) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	6 April 2022

For research use only

Publishing research using A 6841? Please let us know so that we can cite your publication as a reference.

