



Data Sheet

ANTI HUMAN PROLINE-RICH ACIDIC PROTEIN 1 (PRAP1) PEPTIDE (aa 40-51)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1063.1 / AE1063.2
Immunogen:	Synthetic human Proline-rich acidic protein 1 (PRAP1) Peptide (aa 40-51), polyLysine-conjugated (VLFPPVQKPKLL)
Swiss-Prot No:	Q96NZ9
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Proline-rich acidic protein 1 (PRAP1) Peptide (aa 40-51)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:10.000); Western Blot (1:15.000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	12 May 2023

For research use only

Publishing research using AE1063? Please let us know so that we can cite your publication as a reference.

