



Data Sheet

ANTI HUMAN PROLINE-RICH ACIDIC PROTEIN 1 (PRAP1) PEPTIDE (aa 64-73)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|---|
| Catalog no.: | AE1064.1 / AE1064.2 |
| Immunogen: | Synthetic human Proline-rich acidic protein 1 (PRAP1) Peptide (aa 64-73), polyLysine-conjugated (PILPGTKAWM) |
| Swiss-Prot No: | Q96NZ9 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Human Proline-rich acidic protein 1 (PRAP1) Peptide (aa 64-73) |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | ELISA (1:5.000); Western Blot (1:10.000) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.</small> |
| Store at: | 2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 12 May 2023 |

For research use only

Publishing research using AE1064? Please let us know so that we can cite your publication as a reference.

