



Data Sheet

PORCINE RELAXIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 1027.1 / A 1027.2
Immunogen:	Porcine Relaxin
Synonyms:	N/A
Swiss-Prot No:	P01348
Gene Information:	Gene Name: RLN GeneID: 396891
Host:	Rabbit
Matrix:	Serum
Specificity:	Porcine Relaxin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1027? Please let us know so that we can cite your publication as a reference.

