



## Data Sheet

# HUMAN RELAXIN 2 A CHAIN (aa 14-24)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 9221.1 / A 9221.2
<b>Immunogen:</b>	Synthetic human Relaxin 2 A chain (aa 14-24) (GCTKRSLARFC- $\beta$ Ala-C) KLH-conjugated
<b>Swiss-Prot No:</b>	P04090
<b>Gene Information:</b>	Gene Name: RLN2 GeneID: 6019
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Relaxin 2
<b>Contents:</b>	20 $\mu$ l / 100 $\mu$ l (lyophilized) Resuspend in 20 $\mu$ l / 100 $\mu$ l aqua bidest.
<b>Known applications:</b>	ELISA (1:1000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	13 April 2022

### For research use only

**Publishing research using A 9221? Please let us know so that we can cite your publication as a reference.**

