



Data Sheet

HUMAN RELAXIN 1 B CHAIN (aa 1-9)

ANTIBODY, POLYCLONAL

Catalog no.:	A 9230.1 / A 9230.2
Immunogen:	Synthetic human Relaxin 1 B chain (aa 1-9) (VAAKWKDDV- β Ala-C) KLH-conjugated
Swiss-Prot No:	P04808
Gene Information:	Gene Name: RLN1 GenelD: 6013
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Relaxin 1 B chain (aa 1-9) There was no cross reactivity obtained with human Relaxin 1 A chain.
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 April 2022

For research use only

Publishing research using A 9230? Please let us know so that we can cite your publication as a reference.

