



Data Sheet

PORCINE RELAXIN

ANTIBODY, POLYCLONAL

Catalog no.: AS1024.1 / AS1024.2

Immunogen: Porcine Relaxin

Swiss-Prot No: P01348

Gene Information: Gene Name: RLN
GeneID: 396891

Host: Rabbit

Matrix: Serum

Specificity: Porcine Relaxin

Contents: 20 µl / 100 µl (lyophilized)
Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: RIA (1:100 000)¹, immunohistochemistry (paraffin sections, cryosections, 1:2 000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

References: 1. Kotwica G, Dusza L, Ciereszko R, Okrasa S, Schams D (1991). Evidence for relaxin and progesterone synchronous secretion on days 13 to 17 of the oestrous cycle in sows. *Exp Clin Endocrinol* **98**(1): 3-8.

Last updated on: 28 April 2022

For research use only

Publishing research using AS1024? Please let us know so that we can cite your publication as a reference.

