



Data Sheet

HUMAN RELAXIN 1

ANTIBODY, MONOCLONAL

Catalog no.:	A 9215.1 / A 9215.2
Immunogen:	Human Relaxin 1
Swiss-Prot No:	P04808
Gene Information:	Gene Name: RLN1 GenelID: 6013
Host:	Mouse Balb/c
Clone no.:	6F1willi
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Protein G purified, PBS pH 7.4, 0.1% NaN ₃
Specificity:	Human Relaxin 1 There was no cross reactivity obtained with human Relaxin 2.
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA, FACS This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 April 2022

For research use only

Publishing research using A 9215? Please let us know so that we can cite your publication as a reference.

