



## Data Sheet

# HUMAN RELAXIN 2

## ANTIBODY, MONOCLONAL

<b>Catalog no.:</b>	A 9216.1/A 9216.2
<b>Immunogen:</b>	Human Relaxin 2
<b>Swiss-Prot No:</b>	P04090
<b>Gene Information:</b>	Gene Name: RLN2 GeneID: 6019
<b>Host:</b>	Mouse Balb/c
<b>Clone no.:</b>	2F1fred
<b>Isotype:</b>	IgG <sub>1</sub>
<b>Matrix:</b>	Cell culture supernatant, Protein G purified, PBS pH 7,4, 0,1% NaN <sub>3</sub>
<b>Specificity:</b>	Human Relaxin 2 There was no cross reactivity obtained with human Relaxin 1.
<b>Contents:</b>	10 µg / 100µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA, FACS This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	23 September 2024

### For research use only

**Publishing research using A 9216? Please let us know so that we can cite your publication as a reference.**

