

Data Sheet

HUMAN RELAXIN-2

ANTIBODY, POLYCLONAL

Catalog no.:	A 9231.3 / A 9231.4
Immunogen:	Synthetic human Relaxin-2
Swiss-Prot No:	P04090
Gene Information:	Gene Name: RLN2 GenelD: 6019
Host:	Rabbit
Matrix :	IgG, Protein A purified from serum; PBS pH 7.4, 20 mM D-(+)-trehalose dihydrate, 40 mM D-(-)-mannitol
Specificity:	Human Relaxin-2
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA (2.5 µg/ml) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 March 2025

For research use only

Publishing research using A 9231? Please let us know so that we can cite your publication as a reference.

