



Data Sheet

ANTI HUMAN PRO-RELAXIN 1/2

ANTIBODY, MONOCLONAL

Catalog no.:	A 9214.1 / A 9214.2
Immunogen:	Human pro-Relaxin 1/2
Host:	Mouse Balb/c
Clone no.:	kurd
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Protein A purified, PBS pH 7.4, 0.09% NaN ₃
Specificity:	Human pro-Relaxin 1/2
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA, FACS This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 April 2022

For research use only

Publishing research using A 9214? Please let us know so that we can cite your publication as a reference.

