

Data Sheet

HUMAN RELAXIN 2 A CHAIN (aa 14-24)

SYNTHETIC

Catalog no.:	A 9221AG.1 / A 9221AG.2
Name:	Human Relaxin 2 A chain (aa 14-24)
Sequence:	GCTKRSLARFC- β Ala-C
Swiss-Prot No:	P04090
Gene Information:	Gene Name: RLN2 GeneID: 6019
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9221AG? Please let us know so that we can cite your publication as a reference.

