

Data Sheet

HUMAN RELAXIN-2

SYNTHETIC

Catalog no.:	A 9210AG.1 / A 9210AG.2
Name:	Human Relaxin-2, trifluoroacetate salt <u>Fully native sequence:</u> Relaxin A chain (QLYSALANKCCHVGCTKRSLARFC), Relaxin B chain (DSWMEEVIKLCGRELVRAQIAICGMSTWS) <i>Modified residues Q1': L-pyroglutamic acid</i>
MW:	5962.7 Da
Swiss-Prot No:	P04090
Gene Information:	Gene Name: RLN2 GeneID: 6019
Source:	Synthetic
Purity:	> 98% (HPLC)
Contents:	20 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	06 May 2022

For research use only

Publishing research using A 9210AG? Please let us know so that we can cite your publication as a reference.