

## Data Sheet

# HUMAN RELAXIN 2 A CHAIN (aa 3-10)

## SYNTHETIC

<b>Catalog no.:</b>	A 9220AG.1 / A 9220AG.2
<b>Name:</b>	Human Relaxin 2 A chain (aa 3-10)
<b>Sequence:</b>	YSALANKC- $\beta$ Ala-C
<b>Swiss-Prot No:</b>	P04090
<b>Gene Information:</b>	Gene Name: RLN2 GeneID: 6019
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	11 May 2022

### For research use only

**Publishing research using A 9220AG? Please let us know so that we can cite your publication as a reference.**

