



## Data Sheet

# HUMAN RESISTIN (aa 30-49)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 8027.1 / A 8027.2
<b>Immunogen:</b>	Synthetic human Resistin (aa 30-49) (IGLECQSVTSRGLATCPRG)
<b>Synonyms:</b>	Cysteine-rich secreted protein FIZZ3, Adipose tissue-specific secretory factor (ADSF), C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2
<b>Swiss-Prot No:</b>	Q9HD89
<b>Gene Information:</b>	Gene Name: RETN, FIZZ3, HXCP1, RSTN GeneID: 56729
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Resistin peptide (aa 30-49); human Resistin
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	08 April 2022

### For research use only

**Publishing research using A 8027? Please let us know so that we can cite your publication as a reference.**

