



Data Sheet

HUMAN RESISTIN (aa 64-88)

ANTIBODY, POLYCLONAL

Catalog no.:	A 8026.1 / A 8026.2
Immunogen:	Synthetic human Resistin (aa 64-88) bTG- conjugated (DVRAETTCHCQCAGMDWTGARCCRV)
Synonyms:	Adipose tissue-specific secretory factor, ADSF, C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2, Cysteine-rich secreted protein FIZZ3
Swiss-Prot No:	Q9HD89
Gene Information:	Gene Name: RETN, FIZZ3, HXCP1, RSTN GenelD: 56729
Host:	Rabbit
Matrix:	PBS
Specificity:	Human Resistin
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	08 April 2022

For research use only

Publishing research using A 8026? Please let us know so that we can cite your publication as a reference.

