



Data Sheet

HUMAN RESISTIN (aa 64-88)

ANTIBODY, POLYCLONAL

Catalog no.:	A 8028.1 / A 8028.2
Immunogen:	Synthetic human Resistin (aa 64-88) bTG- conjugated (DVRAETTCHCQCAGMDWTGARCCRV)
Synonyms:	Adipose tissue-specific secretory factor, ADSF, C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2, Cysteine-rich secreted protein FIZZ3
Swiss-Prot No:	Q9HD89
Gene Information:	Gene Name: RETN, FIZZ3, HXCP1, RSTN GeneID: 56729
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Resistin peptide (aa 64-88); human Resistin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 April 2022

For research use only

Publishing research using A 8028? Please let us know so that we can cite your publication as a reference.

