

## Data Sheet

# HUMAN RESISTIN (aa 30-49)

## SYNTHETIC

<b>Catalog no.:</b>	A 8027AG.1 / A 8027AG.2
<b>Name:</b>	Human Resistin (aa 30-49)
<b>Synonyms:</b>	Cysteine-rich secreted protein FIZZ3, Adipose tissue-specific secretory factor (ADSF), C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2
<b>Sequence:</b>	IGLECSVTSRGLATCPRG
<b>MW:</b>	
<b>Swiss-Prot No:</b>	Q9HD89
<b>Gene Information:</b>	Gene Name: RETN, FIZZ3, HXCP1, RSTN Gene ID: 56729
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	6 May 2022

### For research use only

**Publishing research using A 8027AG? Please let us know so that we can cite your publication as a reference.**

