



## Data Sheet

# HUMAN BIG RESISTIN (aa 51-108)

## SYNTHETIC

|                          |  |
|--------------------------|--|
| <b>Catalog no.:</b>      | A 8030AG.1 / A 8030AG.2  |
| <b>Name:</b>             | Human Big Resistin (aa 51-108)   |
| <b>Synonyms:</b>         | Adipose tissue-specific secretory factor, ADSF, C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2, Cysteine-rich secreted protein FIZZ3 |
| <b>Swiss-Prot No:</b>    | Q9HD89   |
| <b>Gene Information:</b> | Gene Name: RETN, FIZZ3, HXCP1, RSTN<br>Gene ID: 56729  |
| <b>Source:</b>           | Synthetic  |
| <b>Purity:</b>           | > 95% (HPLC)   |
| <b>Contents:</b>         | 10 µg / 100 µg (lyophilized)<br>Resuspend in aqua bidest for required concentration.   |
| <b>Store at:</b>         | - 20 °C<br>Repeated thawing and freezing must be avoided   |
| <b>Last updated on:</b>  | 6 May 2022   |

### For research use only

**Publishing research using A 8030AG? Please let us know so that we can cite your publication as a reference.**

